

Alpha Image Reconstruction (AIR) Optimization and Image Quality Assessment

Sergej Lebedev, Stefan Sawall, Stefan Kuchenbecker,
Sebastian Faby, Michael Knaup, and Marc Kachelrieß

German Cancer Research Center (DKFZ), Heidelberg, Germany



DEUTSCHES
KREBSFORSCHUNGSZENTRUM
IN DER HELMHOLTZ-GEMEINSCHAFT

Aims

- Increase convergence speed of the AIR¹ algorithm.
- Demonstrate that conventional image quality metrics can be applied to the AIR images.

¹Hofmann, Sawall, Knaup, and Kachelrieß, "Alpha Image Reconstruction (AIR): A New Iterative CT Image Reconstruction Approach Using Voxel-Wise Alpha Blending", *Med. Phys.* 41(6), p. 061914 (14 pages), 2014

AIR¹

- AIR = Iterative algorithm that combines multiple, voxel-wise weighted basis images into one new image f_{AIR}

$$f_{\text{AIR}} = \sum_{b=1}^B \alpha_b \circ f_b$$

Contains weighting coefficients for every voxel

α_b = weighting images

f_b = basis images

$\alpha \circ f$ = Hadamard product


B = number of basis images

- Idea: each basis image has different characteristics such as low noise or high spatial resolution → all characteristics combined in one new image.

¹Hofmann, Sawall, Knaup, and Kachelrieß, "Alpha Image Reconstruction (AIR): A New Iterative CT Image Reconstruction Approach Using Voxel-Wise Alpha Blending", *Med. Phys.* 41(6), p. 061914 (14 pages), 2014

AIR

- Basic application: combination of a sharp and smooth FBP.



AIR

- AIR minimizes a cost function:

$$C(\boldsymbol{\alpha}) = \underbrace{\|\mathbf{X} \left(\sum_{b=1}^B \boldsymbol{\alpha}_b \circ \mathbf{f}_b \right) - \mathbf{p}\|_{\mathbf{W}}^2}_{\text{convex}} + \underbrace{\beta \sum_{b=1}^B TV(\boldsymbol{\alpha}_b)}_{\text{TV regularization, convex}} + \underbrace{\gamma \sum_{b=1}^B \|\boldsymbol{\alpha}_b - \mathbf{d}_b\|_2^2}_{\text{Penalizes deviations from default images } \mathbf{d}_b, \text{ strictly convex.}}$$


- Gradient descent approach:

$$\boldsymbol{\alpha}_b^{\nu+1} = \boldsymbol{\alpha}_b^{\nu} - \lambda \nabla_{\boldsymbol{\alpha}_b} C(\boldsymbol{\alpha}_b^{\nu}) \quad \nabla C(\boldsymbol{\alpha}_b) = \mathbf{f}_b \circ \left(\mathbf{X}^T \mathbf{W} \left(\mathbf{X} \left(\sum_{b=1}^B \boldsymbol{\alpha}_b \circ \mathbf{f}_b \right) - \mathbf{p} \right) \right) + \nabla_{\boldsymbol{\alpha}_b} \left(\beta \sum_{b=1}^B TV(\boldsymbol{\alpha}_b) + \gamma \sum_{b=1}^B \|\boldsymbol{\alpha}_b - \mathbf{d}_b\|_2^2 \right)$$

- $\beta, \gamma = 0.01$

Improved AIR

- Improved AIR separates optimization into two steps:



- $\beta, \gamma = 0.01$ (different weighting between the penalty terms is possible)

Improved AIR

- Improved AIR algorithm uses a method similar to the improved total variation (iT^V¹) method.
 - A new condition is introduced with which a new τ is determined automatically in each iteration.
- Condition for the rawdata fidelity after each iteration:

$$\underbrace{\|\mathbf{X}\mathbf{f}^n - \mathbf{p}\|_2^2}_{\text{Rawdata fidelity of the image from current iteration}} = (1 - \omega) \underbrace{\epsilon_{SART}^n}_{\text{SART rawdata fidelity from current iteration}} + \omega \underbrace{\epsilon^{n-1}}_{\text{Rawdata fidelity: image from previous iteration}}, \quad \omega \in [0, 1]$$

with

$$\mathbf{f}^n = \sum_{b=1}^B ((1 - \tau)\alpha_{b,SART} + \tau\alpha_{b,Reg}) \circ \mathbf{f}_b$$

¹Ritschl, Bergner, Fleischmann and Kachelrieß, "Improved total variation-based CT image reconstruction applied to clinical data", *Phys. Med. Biol.* 56(6), pp. 1554-1561, 2011

Improved AIR

- Linear combination of the SART and the regularized image:

$$\alpha = (1 - \tau)\alpha_{SART} + \tau\alpha_{Reg}, \quad \tau \in [0, 1]$$

$$\tau_{1,2} = A \pm \sqrt{A^2 + B}$$

$$A = \frac{\mathbf{X}f_{SART}^n(\mathbf{X}f_{SART}^n - \mathbf{X}f_{Reg}^n - \mathbf{p}) + \mathbf{X}f_{Reg}^n\mathbf{p}}{(\mathbf{X}f_{SART}^n - \mathbf{X}f_{Reg}^n)^2}$$



$$B = \omega(\epsilon^{n-1} - \epsilon_{SART}^n)$$

f_{SART}/f_{Reg} = AIR CT images calculated using only $\alpha_{SART}/\alpha_{Reg}$

- Parameter ω chosen manually/globally.
 - Sets rawdata fidelity value for each iteration
 - Controls the strength of the regularization
- $\omega=0.2$ was chosen.

Performance


- High quality images can be acquired after a couple 1000 iterations of the gradient descent implementation and after 200-300 iterations of the improved algorithm.



Convergence plots of the rawdata fidelity


AIR Applied to Clinical Data

FBP(B10f-Kernel)
 $\sigma=32$ HU



C/W = 0/500 HU


FBP(B50f-Kernel)
 $\sigma=111$ HU




C/W = 0/500 HU

AIR Applied to Clinical Data


FBP (B50f)
 $\sigma=111$ HU



α_1




AIR
 $\sigma=54$ HU



FBP (B10f)
 $\sigma=32$ HU



α_2



C/W = 0/500 HU

C/W = 0.5/1

C/W = 0/500 HU

$\omega=0.2$

AIR Applied to Clinical Data

FBP (B50f-Kernel)




AIR




C/W = 0/500 HU

$\omega=0.2$



AIR Applied to Clinical Data


Smooth FBP
(Filtered with a Gaussian filter)
 $\sigma=6$ HU



C/W = 0/500 HU

C/W = 0.5/1

FBP(B50f-Kernel)
 $\sigma=79$ HU



C/W = 0/500 HU


Noise levels in soft tissue (yellow circle): FBP(B50f): $\sigma=79$ HU, smooth FBP: $\sigma=6$ HU

AIR Applied to Clinical Data

FBP (B50f)
 $\sigma=79$ HU



α_1





AIR
 $\sigma=43$ HU

Smooth FBP
 $\sigma=6$ HU



α_2



C/W = 0/500 HU

C/W = 0.5/1

C/W = 0/500 HU


Noise levels in soft tissue (yellow circle): FBP(B50f): $\sigma=79$ HU, smooth FBP: $\sigma=6$ HU

AIR: $\sigma=43$ HU


$\omega=0.2$

AIR Applied to Clinical Data

FBP (B50f-Kernel)



AIR



C/W = 0/500 HU

$\omega=0.2$

Modulation Transfer Function


- If an MTF or another image quality metric is defined for the basis images it can be estimated for every voxel of the AIR image.
- MTF of a B50f/B10f-Kernel was measured at a Definition Flash Scanner.

$$MTF(j, \rho) = \sum_n^B \alpha_b^j MTF_b(\rho)$$

α_b^j = voxel j of weighting image b

B = number of basis images


MTF_b = MTF of the basis image b



Modulation Transfer Functions for the Basis Images

Modulation Transfer Function

- The MTF is computed for each voxel.
- The “10%-value“ of the MTF ρ_{10} for each voxel is displayed as a map.



Conclusion

- **Optimized AIR algorithm improves performance by a factor of about 5-10.**
- **Noise can be significantly reduced while spatial resolution at edges is mostly maintained.**
- **Predictions for image quality metrics based on the basis images are possible.**

Thank You!



The 4th International Conference on
Image Formation in X-Ray Computed Tomography

July 18 – July 22, 2016, Bamberg, Germany
www.ct-meeting.org



© Archiv des Bamberg Tourismus & Kongress Service

Conference Chair

Marc Kachelrieß, German Cancer Research Center (DKFZ), Heidelberg, Germany

This presentation will soon be available at www.dkfz.de/ct.

Parts of the reconstruction software RayConStruct IR were provided by
RayConStruct[®] GmbH, Nürnberg, Germany.