

PROJECT WAS DISCUSSED WITH MCTP-STAFF

metabolomics-projects@uni.heidelberg.de

ONLINE REQUEST SUBMITTED (PER INDIVIDUAL ANALYSIS) hmls.corefacilities.org/

## SAMPLES PREPARED ACCORDING TO USER GUIDELINES

- Initials + consecutive numbering as label, e.g. MM\_01
  - No adhesives only written labels
  - Sample sets >15 samples, please randomize sample order
- For correct sample weights/volumes, please check user guidelines
  - Solids: Weighed and pulverized in 1.5/2ml tubes
  - Solids whole: weighed sample in 2ml tubes
  - Liquids: Not overfilled in 1.5/2ml tubes

## SAMPLE HANDOVER

- Transfer on dry ice
- Delivery in plastic or paper box

Submission times: MON 10:30 – 11:30, TUE 9:00 – 11:00, WED 13:30 – 15:00, THU 9:00 – 11:00; other times upon prior consultation via metabolomics-samples@uni.heidelberg.de

